

LABORATORY CHEMICAL FUME HOOD INSPECTION	DATE OF PREVIOUS INSPECTION _____	DATE _____
	THIS INSPECTION PERFORMED BY (Name) _____	
LOCATION OF HOOD _____	TYPE OF HOOD <input type="checkbox"/> Standard <input type="checkbox"/> Auxiliary Air supply <input type="checkbox"/> Other (specify) _____	
GENERAL TOXICITY RATING OF MATERIAL USED IN HOOD <input type="checkbox"/> Low (STEL > 1,000 PPM) <input type="checkbox"/> Medium <input type="checkbox"/> High (STEL < 10 PPM)	CROSS SECTIONAL AREA AT FACE Height: _____ feet x Width: _____ feet = _____ feet ²	


AIR VELOCITY READINGS
(Readings are to be taken in the center at each of the prescribed frontal grids.)

FPM	FPM	FPM
FPM	FPM	FPM
FPM	FPM	FPM

Exhaust on, sash fully raised.
(Exhaust flow value equal to zero CFPM)

$$\frac{\text{FPM} + \text{FPM} + \text{FPM} + \dots + \text{FPM}}{9} = \text{FPM average.}$$

Average value _____ FPM.



FPM	FPM	FPM
FPM	FPM	FPM

Exhaust on, sash raised 18 inches.
(Readings may not vary more than ± 10 FPM from average value.)

Average value _____ FPM.
(Value should be 80- 100 FPM.)

Exhaust flow value _____ CFPM.

FPM	FPM	FPM
-----	-----	-----

Exhaust on, sash 6 inches above work surface.
(Readings shall be at least 2 but not more than 3 times the face velocity when sash was fully raised)

Average value _____ FPM.

Exhaust flow value _____ CFPM

EXHAUST READING WITH SASH CLOSED

Exhaust flow value _____ CFPM.

TITANIUM TETRACHLORIDE INDICATION OF FLOW PATTERNS AT HOOD FACE.

- ☐ Satisfactory flow patterns evident.
☐ Unsatisfactory (describe): _____

ONE-MINUTE SMOKE BOMB DISCHARGE

- ☐ Effective smoke removal with sash fully raised.
☐ Effective smoke removal with sash 6 inches above work surface.
☐ Effective smoke removal with sash closed.
☐ If unsatisfactory, describe: _____

APPROVAL

- ☐ This hood is found to be acceptable for use with materials of the general toxicity rating as specified above.
☐ This hood has been found UNACCEPTABLE.

SIGNATURE _____

DATE _____

HOW TO MEASURE FUME HOOD FACE VELOCITY: STANDARDIZED FOR REE INCLUDING THE SHEMA AND ENGINEERING COMMUNITIES

Note: The sole purpose of this document is for it to serve as a user-friendly and quick-reference guide on how to measure the airflow of a chemical fume hood. It is not meant to replace information currently contained in ARS Manual 230. Always refer to M230, the American National Standards for Laboratory Ventilation (ANSI/AIHA Z9.5-1992), the ASHRAE Guideline: Method of Testing Performance of Laboratory Fume Hoods (ANSI/ASHRAE 110-1995), the ASHRAE 1995 HVAC Applications Handbook, Fire Protection for Laboratories Using Chemicals (NFPA 45), Laboratory Fume Hoods Recommended Practices (SEFA 1.2-1996) and 29 CFR 1910.1000 for detailed specific information.

1. To measure chemical fume hood face velocities use a calibrated Thermal Anemometer, Velometer, or Swinging Vane Anemometer.
2. It is recommended that you mount the measuring equipment to a sturdy, freestanding device such as a ring stand to allow for steady & accurate airflow measurement (this helps to negate the possibility of generating invalid data due to unsteady hand-held equipment).
3. Close all doors and windows and minimize foot traffic in the work area.
4. While the fume hood is running measure the Cross-Draft Velocity (airflow perpendicular to fumehood intake) six inches in front of the mid-point of the hood. This reading should not exceed 20 FPM.
5. Fully open the sash.
6. Visually check that air is flowing into the fumehood and baffles using a smoke-generating device (titanium tetrachloride smoke sticks, sulfuric acid mist tubes, or dry ice). No smoke should escape from the hood.
7. Divide the wide opening into nine equally spaced sectors (3 columns/3 rows). Refer to USDA S&E Form 283.
8. Move the measuring device to the center of each of the 9 sectors and record the face velocity.
9. No face velocity measurement on the grid should be less than 60 FPM.
10. Average the 9 face velocity measurements—this average should be 80-120 FPM.
11. Lower the sash to 18 inches above the work surface.
12. Divide the opening into six equally spaced sectors (3 columns/2 rows). Refer to Form 283.
13. Move the measuring device to the center of the 6 sectors and record the face velocity.
14. Average the 6 face velocity measurements—this average should be 100-150 FPM.
15. Lower the sash to 6 inches above work surface.
16. Divide this opening into three equally spaced sectors (3 columns/1 rows). Refer to Form 283.
17. Move the measuring device to the center of each of the 3 sectors and record the face velocity.
18. Average the 3 face velocity measurements—this average should be 160-300 FPM.
19. To minimize turbulence and reverse flow the Average Face Velocity should not exceed 150 FPM at the working sash level (18 inches).
20. The Cross-Draft Velocity should not exceed 20% of the Average Face Velocity.
21. Adjust the fumehood baffles or HVAC system to ensure items 19 & 20.
22. Record the measurements and label the fumehood with the measurement data (Average Face Velocity, Sash Height, Date, & Inspector).

The following three conditions must be met in order for a hood to pass:

1. The average face velocity at a sash opening of 18 inches must be a minimum of 100 FPM.
2. The average face velocity at a sash opening of 6 inches cannot exceed 300 FPM.
3. At any sash height smoke cannot escape from the hood back into the room.